

**Valoroux 3Ply Disposable Face Masks**

<b>Test</b>	<b>ASTM F2100</b>
	<b>Equivalent to Level 1</b>
<b>Bacterial Filtration Efficiency, %</b>	≥95
<b>Differential Pressure mm H<sub>2</sub>O/cm<sup>2</sup> (Pa/cm<sup>2</sup>)</b>	<5.0
<b>Submicron particulate filtration efficiency at 0.1 micron, %</b>	≥95
<b>Splash Resistance/Synthetic Blood Resistance, mmHg</b>	80
<b>Flame Spread</b>	Class 1
<b>Microbial Cleanliness (cfu/g)</b>	N/A

<b>Face Fit</b>	Loose-fitting
<b>Intended use</b>	For use as a source control by the public to help prevent the spread of infection or illness
<b>Filtration Efficiency</b>	Degree of protection equivalent to ASTM Level 1
<b>Air Leakage</b>	At the edges/sides of mask when user inhales
<b>Tests and Approvals</b>	Testing done at Health Canada approved Lab, Medical Device licence

**How to wear a mask**

A mask can provide the protection if it is worn properly.

This means there should be no gap along the edges of the mask between the face and the mask and that the mask should cover from the bridge of the nose to the jawline/chin. It is vital that the mask covers not only the mouth but also the nose. To guarantee adequate coverage of the nose, masks should not be worn below or at the tip of the nose.

**Additional information**

Model Number	VA-BHT-001C
Color	Blue
ASTM level	Equivalent to Level 1
Packaging Detail	50/box